

Speed-dating with Praat ... exploring some of its basic functions!

Youri Maryn, PhD

PO = window "Praat Objects"

PE = window "Praat Edit"

PP = window "Praat Picture"

2 – Copy, listen to, downsample and re-listen to a sound object

2.1	PO	Copy...
2.2	Copy object	Type another name in the text line OK
2.3	PO	Play
2.4	PO	Convert -
2.5	PO	Resample...
2.6	Sound: Resample	New sampling frequency (Hz): 26000 Precision (samples): 50 OK
2.7	PO	Play
2.8	Sound: Resample	New sampling frequency (Hz): 10000 Precision (samples): 50 OK
2.9	PO	Play
2.10	Sound: Resample	New sampling frequency (Hz): 5000 Precision (samples): 50 OK
2.11	PO	Play